90-56725-03

August 2014

## RESTYLANE FINE LINES-L™ INSTRUCTIONS FOR USE

**Formula**  
Contents:  
Hyaluronic acid, stabilized 20 mg/ml  
Lidocaine hydrochloride 0.3%  
Phys. sodium chloride solution, pH 7.4 s.

**Description**  
*Restylane Fine Lines-L* is a clear, transparent and viscous gel supplied in a glass syringe together with a 30 G or 29 G Thin Wall (TW) needle. The product is for single use only. *Restylane Fine Lines-L* is a unique form of non-animal, stabilized hyaluronic acid (NASHA™). Hyaluronic acid is a natural polysaccharide which occurs as an important structural element in the skin and in subcutaneous and connective tissues as well as in the synovial tissue and fluid. Hyaluronic acid belongs to a group of very few substances which are identical in all living organisms.

**Mode of action**  
*Restylane Fine Lines-L* acts by adding volume to the tissue, thereby restoring the skin contours to the desired level of correction. *Restylane Fine Lines-L* is naturally integrated into the tissue and will in time undergo isovolemic degradation.

*Restylane Fine Lines-L* has been enhanced with the addition of lidocaine in order to reduce patient discomfort during treatment.

**Indication and usage**  
*Restylane Fine Lines-L* is intended to be used for facial tissue augmentation. It is recommended to use the product for the correction of thin superficial lines, for example periorbital lines, perioral lines, forehead lines and smile lines. It should be injected in the superficial dermal layer of the facial skin. Deeper injections into the subcutaneous fatty tissue or supraperiosteal layer are appropriate for areas with adequate soft tissue support and soft tissue cover such as midface and jaw line. With cutaneous contour deformities the best results are obtained if the defect can be manually stretched to the point where it is eliminated. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. Markedly indurated defects may be difficult to correct. For the correction of cutaneous contour deformities such as wrinkles it is recommended that *Restylane®* be used. *Perlane®* can be used for shaping the contours of the face and for the correction of folds using the co-packed needles or a small gauge blunt cannula. The blunt cannula is appropriate for areas with adequate soft tissue support and soft tissue cover such as the midface and jaw line. For lip enhancement both *Restylane* and *Perlane* can be used. Also combinations of *Restylane*, *Perlane* and *Restylane Fine Lines™* can be used. Please consult the *Restylane* and *Perlane* Instructions for Use for more information.

**Contraindications**  
Do not use *Restylane Fine Lines-L* on patients with known hypersensitivity to lidocaine or amide type local anesthetics.

**Warning**  
*Restylane Fine Lines-L* is only intended for use as an intradermal implant. Do not inject intravascularly. Do not mix with other products.

**Precautions**  
Normal precautions associated with intradermal injections must be observed. Like any such procedure, the implantation of *Restylane Fine Lines-L* is associated with an inherent risk of infection. *Restylane Fine Lines-L* should not be used in or near anatomic sites where there is active skin disease, inflammation or related conditions. Do not use *Restylane Fine Lines-L* together with any other injectable implant, except for *Restylane* and *Perlane* products.

*Restylane Fine Lines-L* should not be injected into an area where a permanent implant has been placed.

*Restylane Fine Lines-L* should not be used for patients with unattainable expectations. The patient should be informed that he or she should not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold at least until the initial swelling and redness have resolved.

If laser treatment, chemical peeling or any other procedure based on active dermal response is considered after treatment with *Restylane Fine Lines-L* there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if *Restylane Fine Lines-L* is administered before the skin has healed completely after such a procedure.

There is a potential risk with the procedure that the material could be inadvertently injected into dermal blood vessels, which could lead to vascular occlusion in an end-artery with corresponding consequences. *Restylane Fine Lines-L* has not been tested in pregnant or lactating women or in children.

**Anticipated side-effects**  
After the injection of *Restylane Fine Lines-L*, some common injection-related reactions might occur. These reactions include erythema, swelling, pain, itching, discoloration or tenderness at the implant site. Typically resolution is spontaneous within a few days after injection into the skin. Additionally, temporary palpable lumpiness has been noted after the use of *Restylane Fine Lines-L* in some patients. The composition of *Restylane Fine Lines-L* is set to balance a normal tissue pressure. However, because the tissue pressure is sometimes disturbed to a higher value such as during oedema or to a lower value such as during dehydration, a small but significant change (swelling or shrinkage) may occur.

**Adverse events**  
Reactions thought to be of a hypersensitivity nature have been reported in about 1 in every 10,000 treatments. These have consisted of swelling and induration at the implant site, sometimes with oedema in the surrounding tissues. Erythema, tenderness and rarely acneiform papules may also occur. The reactions have started either shortly after injection or after a delay of 2-4 weeks and have generally been described as mild to moderate and self-limiting, with an average duration of 2 weeks. In pronounced cases a short course of oral corticosteroids may prove effective. In addition, rare cases (less than 1 in 15,000 treatments) of granuloma formation, superficial necrosis and urticaria have been reported. Patients who have experienced this type of reaction should not be retreated with a *Restylane* product.

Adverse events must be reported to Galderma Canada Inc.

**Interactions**  
Treatment with *Restylane Fine Lines-L* in combination with other drugs and devices has not been tested.

**Assembly of 30 G needle to syringe**  
For safe use of *Restylane Fine Lines-L* it is important that the needle is properly assembled.


**See pictures A and B.**

A. Unscrew the tip cap of the syringe carefully.


B.1. Take a loose grip on the needle shield and mount the needle on the Luer-Lok® by screwing until you feel some counterpressure.

B.2. Take a firmer grip, **press and turn** a further 90° (a quarter of a turn).

B.3. Pull off the needle shield.

**Assembly of 29 G (TW) needle to syringe**  
Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the Luer-Lok adapter. Grasp the needle shield with the other hand. To facilitate proper assembly, both push and rotate firmly.

**Disposal and administration**  
Before the treatment, the patient's suitability for the treatment and the need for pain relief should be assessed. Normally, no anesthesia is necessary when treating thin superficial lines. The patient should be informed about the indications, expected result, contraindications, precautions, warnings and potential adverse events. The treatment site should be cleaned with a suitable antiseptic solution. *Restylane Fine Lines-L* is administered using a thin gauge needle (30 G or 29 G (TW)) by injecting the material into the superficial layer of the skin. An injection too superficial may give blanching effects and bumps on the treatment site. If blanching is observed, i.e. the overlying skin turns a whitish colour, the injection should be stopped at once and the area massaged until it returns to a normal colour. Before injecting, press the rod carefully until a small droplet is visible at the tip of the needle.

The injection technique with regard to the depth of injection and the administered quantity may vary. The linear threading technique can be used to carefully lift up the thin superficial line, but some physicians prefer a series of punctual injections or a combination of the two. During injection it is recommended that the eye of the needle should face upwards. Inject *Restylane Fine Lines-L* while pulling the needle slowly backwards. Injection should stop just before the needle is pulled out from the skin to prevent material from leaking out from the injection site. Please consult Galderma Canada Inc. for details. Defects should be fully corrected, but not overcorrected, at each treatment session. If the skin of the patient is very loose, it is recommended that *Restylane Fine Lines-L* is injected on two separate occasions. The correction site should be massaged to conform to the contour of the surrounding tissues. For each treatment site a maximum dosage of 0.8 ml is recommended. If the treated area is swollen directly after the injection, melting ice or a cold pack can be applied on the site for a short period. After the first treatment, additional implantations of *Restylane Fine Lines-L* may be necessary to achieve the desired level of correction. Periodic follow-up injections help sustain the desired degree of correction.

**Note** The correct injection technique is crucial for the final result of the treatment. *Restylane Fine Lines-L* is only intended to be administered by authorized personnel in accordance with local legislation.

**The syringe, the needle and any unused material must be discarded directly after the treatment session.**

**How supplied**  
*Restylane Fine Lines-L* is supplied in a glass syringe with a Luer-Lok fitting. *Restylane Fine Lines-L* is co-packed with a sterilized needle(s) as indicated on the carton, either 30 G x 12" or 29 G (TW) x 12". A patient record label forms part of the syringe label (see picture C). Remove by pulling the flap marked with three small arrows (see picture D). This label is to be attached to patient records to ensure traceability of the product. The contents of the syringe are sterile. The number of units per package and the volume contained in each syringe is as stated on the outer package.

**STERILE NEEDLES**  
• Follow national, local or institutional guidelines for use and disposal of medical sharp devices. Obtain prompt medical attention if injury occurs.  
• To help avoid needle breakage, do not attempt to straighten a bent needle. Discard it and complete the procedure with a replacement needle.  
• Do not re-shield used needles. Recapping by hand is a hazardous practice and should be avoided.  
• Discard unshielded needles in approved sharps collectors.

**Shelf life and storage**  
As indicated on the package. Store up to 25°C (77°F). Protect from sunlight and freezing. Refrigeration not required.

Canadian Patent No. 2,226,488  
**Manufactured by**  
Q-Med AB,  
Seminariégatan 21, SE-752 28 Uppsala, Sweden  
Phone +46(0)18 474 90 00, Fax +46(0)18 474 90 01  
www.q-med.com, e-mail: info.q-med@galderma.com

**Manufactured for**  
Galderma Canada Inc.  
Thornhill, ON L3T 7W3  
1 800 467-2081

**IF THE PACKAGE IS DAMAGED, DO NOT USE**

**Symbols on packaging**

Refer to instructions for use

For single use

Sterile. The contents of the syringe have been sterilized by using moist heat.

Sterile. The needle has been sterilized by using irradiation.

*Restylane*, *Perlane*, *Restylane Fine Lines*, and *Restylane Fine Lines-L* are trademarks of Galderma S.A.

All other trademarks are the property of their respective owners.

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